

BLUES II

Wireless, hands-free and one-touch mobile phone communication!

New and improved Bluetooth neck loop for your mobile phone! Blues II transfers quality sound from the mobile phone to your hearing aids/CI, which has t-coil, so that you can talk hands-free with good amplification without background noise and electronic interference.

A multifunctional control button to answer, end and reject calls makes the neck loop extremely flexible.

Blues II is the world's first Bluetooth neck loop that has the latest, ground-breaking Bluetooth technology. It gives you the opportunity to enjoy all the functions on your mobile phone, e.g. conversations, games, music, film etc, and all with a sensational range of up to 30 meters. It is also easy to connect to your mobile phone or PC.

Blues II transfers all sounds from the computer and enables VoIP communication such as Skype, etc.

Blues II has an integrated Speed-Dial button so that you can easily and quickly call a pre-programmed number.

By using innovative technology, we have made it possible to connect to a headset. This gives you the freedom to choose based on your own needs!

Compatible with hearing aids/CI with "T" or "MT" settings.

Using COMMidt's clip-on microphone makes it possible to wear Blues II under clothing. See illustration to the right.



Scale 1:1



Stylish and the perfect hands-free solution for your mobile phone

- Very easy to use
- Headset output
- Speed-Dial button
- Active noise reduction
- Stepless volume control
- Rechargeable batteries
- Good battery capacity
- Eavesdrop secure
- First class quality sound
- No channels (like FM)
- Discreet and small
- A2DP (Hi-Fi sound)

Art.No: 108-002

Blues II has outstanding battery capacity and we are proud to say that it has the best on the market. With a total of 15 hours talk-time, you can enjoy the neck loop all day without having to charge it. Blues II is easy to charge in just 1.5 hours via the charger, USB cable, Twin Charger Cradle, Audiostreamer or car charger. You can also use the unit while it is charging.



Blues II has a multi-turn neck loop, which means that the sound transferred to the hearing aid(s) / CI is of high quality because the field strength is strong and stable.

This improves speech recognition considerably, because the bandwidth and the dynamics are very good compared to traditional neck loops.

Blues II is designed to fit into a hectic lifestyle. It gives freedom, security and comfort so that you can talk on the phone hands-free, e.g. in the car. Exceptionally user-friendly, crystal clear sound and comfortable design makes Blues II the perfect companion in a busy life.

International charger, clip, USB cable and earbuds are included.

Accessories:

Product #	Product Name	Product Description
957-001	Twin Charger Cradle	Double charger
901-001	External Clip-On Microphone 40 cm	2.5mm jack. Enables to wear Blues II under clothing
907-001	12/24V Car Charger	To charge Blues II in the car
530-001	Earbuds (Headset) with microphone	With 3.5 mm jack
521-002	Headset with mic.	With 3.5 mm jack
930-001	USB Cable	Charge and programme Blues II on the computer
955-002	International charger (mains power)	

Sounds better



Tech. info

- Bluetooth 2.0
- BlueCore 5 Multimedia Engine
- Talk time: up to 15 hours
- Standby time: > 500 hours
- Charging time: 1.5 hours
- Battery type: Rechargeable Li-Ion Polymer 3.7 volt
- Dimension: 64 (L) x 35 (W) x 17 (H) mm (2,5 (L) x 1,4 (W) x 0,7 (H) inch.)
- Weight: approx. 40 grams (approx. 1,4 ounces)
- Headset jack

