

# CIBS PERSONAL LISTENER

## Simple speech amplifier with good sound quality!

Use the new CIBS PERSONAL LISTENER with the Neck Loop or with a headset and the wireless Liberto Bluetooth Microphone to obtain the market's best Personal Listener ready for use. Ideal for hard of hearing people with or without hearing aid(s) / CI

When you place Liberto on your conversation partner or a table, it produces very good S/N ratio because the microphone picks up very little background noise or interference.

The headset is easily connected to the neck loop's headset output at the top of the neck loop.

The sound quality is top class and the range of 10-30 meters also gives you excellent flexibility.

The CIBS PERSONAL LISTENER is equipped with volume control that provides the user the option for easy individual adjustment!

Twin Charger Cradle, clip with string and headset are included



**Wireless  
lip synchronized  
Personal  
Listener**

- Inclusive a combined omni- and directional microphone (Liberto)
- Unique S/N-ratio
- Active noise reduction
- Very easy to use
- Stepless volume control
- Long range
- Unique battery capacity
- Eavesdrop secure
- Simple connection
- No channels (like FM)

### Tech. info, CIBS Neck Loop and Liberto

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Bluetooth 2.0</li> <li>• BlueCore05 Multimedia Engine</li> <li>• Talk time: up to 15 hours</li> <li>• Standby time: &gt; 500 hours</li> <li>• Charging time: 1.5 hours</li> <li>• Rechargeable battery</li> <li>• Battery type: Li-Ion Polymer 3.7 volt</li> <li>• Dimension CIBS Neck Loop:<br/>64 (L) x 35 (W) x 17 (H) mm<br/>(2,5 (L) x 1,4 (W) x 0,7 (H) inch.)</li> </ul> | <ul style="list-style-type: none"> <li>• Dimension Liberto:<br/>73 (L) x 35 (W) x 17 (H) mm<br/>(2,9 (L) x 1,4 (W) x 0,7 (H) inch.)</li> <li>• Weight: approx. 40 grams<br/>(approx. 1,4 ounces)</li> <li>• Headset jack</li> <li>• Simple unit pairing</li> <li>• Headset, hands-free and A2DP<br/>(Hi-Fi) profile</li> </ul> |
|--|--|

Art.No: 108-112